

NRSF2 Design Document Appendix

July 9, 2018

Construct 2 orthogonal q vectors:

$$\begin{aligned}\vec{q}_1 &= R_z(-\frac{2\theta}{2})(010) \\ \vec{q}_2 &= R_z(-\frac{2\theta}{2})(100)\end{aligned}$$

where

- $2\theta_{peak}$ = position of peak (which is in plane)

In order to convert \vec{q} from instrument coordinate to sample coordinate, the rotation will be done on ω , χ and ϕ respectively.

The rotation matrix R_p is defined as

$$R_p = R_x(\phi + 90^\circ) \times R_y(\chi) \times R_z(-\omega)$$

where

- ω = incident angle
- χ = χ rotation about (100) of sample
- ϕ = ϕ rotation about sample normal

The projection is then

$$\begin{aligned}\alpha &= \cos^{-1}(Q)(001)) \\ \beta &= \cos^{-1}(Q)(100)) \\ \alpha &= \cos^{-1}(Q')(001)) \\ \beta &= \cos^{-1}(Q')(100)) \\ \alpha_p &= \alpha - \alpha' \\ \beta_p &= \beta - \beta'\end{aligned}$$